



CRITERION 2 – TEACHING-LEARNING AND EVALUATION

2.3.2 Has provision for the use of Clinical Skills Laboratory and Simulation Based Learning the Institution:

Has Basic Clinical Skills Training Models and Trainers for clinical skills in the relevant disciplines.

Has advanced patient simulators for simulation-based training

Has structured programs for training and assessment of students in Clinical Skills Lab / Simulation centre

Conducts training programs for the faculty in the use of clinical skills lab and simulation methods of teaching-learning


To reduce enormous use of paper and printing the ensure data, sign and a seal by the Competent Authority for all the papers, we have used the Class-3 Digital Signatures where a Registration Authority i.e. Dr. Mahipal Singh, Registrar of our University authenticate the documents and responses claimed in this pdf file.



SHOBHIT UNIVERSITY, Gangoh







Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in

Proof of patient simulators for simulation- based training.



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



CPR Training Mannequin

These manikins simulate blood flow to the brain using LED light feedback. These CPR manikins can help your students learn how to adjust their technique in real-time to deliver CPR correctly. This gives them the confidence and the skills to provide more effective CPR.



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



Nebulizer Training Kit


A nebulizer is a small machine that turns liquid medicine into a mist that can be easily inhaled. You sit with the machine and breathe in the medicine through a connected mouthpiece or facemask.



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



Catheterization Simulator

Trainees can learn and practice proper sequence of catheterization steps. [2] Use of simulation enhances visual and spatial thinking and improves hands-eye coordination



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



IM Injection Simulator

Injection Simulator. Here is the perfect volunteer for training your beginning healthcare students to administer intramuscular injections



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



Bandaging practicing simulator


Ideal for teaching wound management skills, dressing changes, and bandaging techniques



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



IV practice Mannequin Arm


This refers to an artificial **arm** that simulates the physical characteristics of a human **arm**, providing a realistic **training** surface



Tel: +91 7830810052

E-mail: registrargangoh@shobhituniversity.ac.in

U.: www.sug.ac.in



Rt Insertion Stimulator

It is used to aspirate gastric contents in case of hematemesis, intestinal obstruction, post-surgery and gastric acid aspiration for gastric function test